



PUSH BUTTON BLOOD COLLECTION SET

SAFER & EASIER



BEVEL TECHNOLOGY

- Unique bevel design requires lower penetration force*
- More comfortable patient experience

INCREASED FLOW RATE

- Needle Cannula increases flow rate 200% (23G) and 60% (21G)*
- Improve patient comfort: Use 23G instead of 21G without compromising flow rate
- Approved for short-term infusion (<2 hrs)

SAFETY WITH THE PUSH OF A BUTTON

- Single-handed activation retracts the needle instantly from the vein into the needle hub
- Minimizes sharp exposure

*Compared to leading-brand standard cannula

KB Medical Group Inc.


We Focus on a Strong Customer Orientation and Innovative, Top Quality Product Solutions.

www.kbmed.com



SAFER & EASIER

Designed to bring comfort for patients and clinicians



PUSH BUTTON BLOOD COLLECTION AND INFUSION SET

| Code | Description | Color | Box Qty | Case Qty |
|-------------|---|--------|---------|----------|
| SBCNP21K-1 | 21g x 3/4" with 12" tubing and luer adapter | Green | 50 | 200 |
| SBCNP23K-1 | 23g x 3/4" with 12" tubing and luer adapter | Blue | 50 | 200 |
| SBCNP25K-1 | 25g x 3/4" with 12" tubing and luer adapter | Orange | 50 | 200 |
| SBCNP21K-1H | 21g x 3/4" with 12" tubing and luer adapter w/holder attached | Green | 20 | 100 |
| SBCNP23K-1H | 23g x 3/4" with 12" tubing and luer adapter w/holder attached | Blue | 20 | 100 |
| SBCNP25K-1H | 25g x 3/4" with 12" tubing and luer adapter w/holder attached | Orange | 20 | 100 |
| SSVSP21 | 21g x 3/4" with 12" tubing w/o luer adapter | Green | 50 | 200 |
| SSVSP23 | 23g x 3/4" with 12" tubing w/o luer adapter | Blue | 50 | 200 |
| SSVSP25 | 25g x 3/4" with 12" tubing w/o luer adapter | Orange | 50 | 200 |

KB Medical Group Inc.

9650 Research Dr, STE A, Irvine, CA 92618 USA

Web: www.kbmed.com | Tel.: +1 (949) 333-3785

