



Specimen Collection





Safety Blood Collection Needle

KBMED Safety blood collection needle is simple and easy to use. The safety mechanism can be activated immediately after the blood draw, helping to protect you against needlestick injury. It is also offered with a pre-attached holder for added convenience and to help meet OSHA single-use holder standards.

Feature :

- The Safety Shield protects against needle injuries
- Large wing size permits improved positioning during use
- Blood flashback is seen clearly through translucent shield



Type 1 Safety butterfly type



Type 2 Safety needle type



After use



Type 3 Safety pen type



After use

Blood Collection Needle



Type 1 Butterfly type



Type 2 Pen type



Type 3 Hypodermic needle type



Size		Color	Code									
OD(mm)	Gauge		Safety Type		Safety Type with pre-attached holder			Normal Type		Normal Type with pre-attached holder		
			Type 1	Type 3	Type 1	Type 2	Type 3	Type 1	Type 2	Type 1	Type 2	Type 3
0.40	27G	Grey	SBCN27-1	SBCN27-3	SBCN27-1H	SBCN27-2H	SBCN27-3H	BCN27-1	BCN27-2	BCN27-1H	BCN27-2H	BCN27-3H
0.45	26G	Brown	SBCN26-1	SBCN26-3	SBCN26-1H	SBCN26-2H	SBCN26-3H	BCN26-1	BCN26-2	BCN26-1H	BCN26-2H	BCN26-3H
0.50	25G	Orange	SBCN25-1	SBCN25-3	SBCN25-1H	SBCN25-2H	SBCN25-3H	BCN25-1	BCN25-2	BCN25-1H	BCN25-2H	BCN25-3H
0.55	24G	Purple	SBCN24-1	SBCN24-3	SBCN24-1H	SBCN24-2H	SBCN24-3H	BCN24-1	BCN24-2	BCN24-1H	BCN24-2H	BCN24-3H
0.60	23G	Blue	SBCN23-1	SBCN23-3	SBCN23-1H	SBCN23-2H	SBCN23-3H	BCN23-1	BCN23-2	BCN23-1H	BCN23-2H	BCN23-3H
0.70	22G	Black	SBCN22-1	SBCN22-3	SBCN22-1H	SBCN22-2H	SBCN22-3H	BCN22-1	BCN22-2	BCN22-1H	BCN22-2H	BCN22-3H
0.80	21G	Green	SBCN21-1	SBCN21-3	SBCN21-1H	SBCN21-2H	SBCN21-3H	BCN21-1	BCN21-2	BCN21-1H	BCN21-2H	BCN21-3H
0.90	20G	Yellow	SBCN20-1	SBCN20-3	SBCN20-1H	SBCN20-2H	SBCN20-3H	BCN20-1	BCN20-2	BCN20-1H	BCN20-2H	BCN20-3H
1.10	19G	Cream	SBCN19-1	SBCN19-3	SBCN19-1H	SBCN19-2H	SBCN19-3H	BCN19-1	BCN19-2	BCN19-1H	BCN19-2H	BCN19-3H
1.20	18G	Pink	SBCN18-1	SBCN18-3	SBCN18-1H	SBCN18-2H	SBCN18-3H	BCN18-1	BCN18-2	BCN18-1H	BCN18-2H	BCN18-3H



Safety Needle Holder

KBMED Safety Needle Holder is applied in the clinical venous sampling which is achieved by blood collection needles to secure the operators



SNH01




SNH02

Needle Holder

NH01

· KBMED single-use holders are compatible with all Vacutainer tubes and blood collection needles



NH02

· The KBMED Blood Transfer holder is used for needleless specimen transfer from a syringe to an evacuated tube or blood culture bottle.



NH03

· The KBMED Luer Lock Access holder is compatible with female luer connections or IV ports designed for luer access.

Blood Collection Tube

We provide a full offering of blood collection tubes in both glass and plastic.

Type	Code
Glass blood collection tube	BCT01
Plastic blood collection tube	BCT02
Micro collection tube	BCT03
Non-vacuum blood collection tube	BCT04
Capillary blood collection tube	BCT05



Category	Item	Additive	Tube Material	Closure Cap					Test Item
				Type C	Tube Size (mm)	Type H	Tube Size (mm)	Color	
Serum Blood Collection Tube	Plain Tube	Plain	Glass/Plastic		13*75		16*100	Red	Clinical Biochemistry, Immunology and Serology Test
	Pro-coagulation Tube	Clot Activator	Glass/Plastic		13*100				
	Gel & Clot Activator Tube	Gel & Coagulant	Glass/Plastic		13*75 13*100 16*100		16*100	Yellow	
Whole Blood Collection Tube	EDTA Tube	Sprayed K2EDTA/ Sprayed K3EDTA	Glass/Plastic		13*75 13*100 16*100		16*100	Purple	Haematology Test (Blood Routine Examination)
	ESR Tube	3.8% Sodium Citrate Buffer(0.129mol/L)	Glass/Plastic Glass		13*75		8*120		Erythrocyte Sedimentation Rate Test
Plasma Blood Collection Tube	Coagulation Tube	3.2% Sodium Citrate Buffer(0.109mol/L)	Glass/Plastic		13*75 13*100		10.25*64	Blue	Coagulation Function Test
	Heparin Tube	Sodium Heparin/ Lithium Heparin	Glass/Plastic		13*75 13*100 16*100			Green	Clinical Chemistry for Emergency Treatment, Blood Rheology Test
	Gel & Heparin Tube	Gel & Sodium Heparin/ Gel & Lithium Heparin	Glass/Plastic Glass/Plastic		13*75 13*100 16*100		16*100	Green	
	Glucose Tube	Sodium Fluoride & Sodium Heparin/ Sodium Fluoride & EDTA/ Sodium Fluoride & Potassium Oxalate	Glass/Plastic		13*75 13*100			Grey	Glucose and Lactate Test
	EDTA & Gel Tube	Gel & Sprayed K2EDTA/ Gel & Sprayed K3EDTA	Glass/Plastic		13*75 13*100 16*100		16*100	Purple	Molecular Biology Test (Such as PCR)

Vacuum Urine Tube



Tube Size (mm)	Additive	Code	Additive	Code
16*100	Boric Acid	UCTB16100	No	UCT16100
13*100	Boric Acid	UCTB13100	No	UCT13100
13*75	Boric Acid	UCTB13075	No	UCT13075



Urine Transfer Device

Code: UCTD120, Size: 100-120ml

Urine Cup



95KPA Compliant

Code: UC90-B2



Regular type

Code: UC120-B1

Temperature label is available




Size	Code	
	Regular Type	95KPA Compliant
20ml	UC20-B1	UC20-B2
40ml	UC40-B1	UC40-B2
60ml	UC60-B1	UC60-B2
90ml	UC90-B1	UC90-B2
120ml	UC120-B1	UC120-B2





Accessories



Tourniquet

Sizes: 1"x18"x0.025"/1"x18"x0.029"/etc.
Latex free, Blue/Orange/Green/Yellow/Pink



Adhesive bandage

Sizes: 72x19mm/76x25mm/Dia. 25mm/etc.
Sizes can be customized



Alcohol Prep Pads

Sizes: 3x6cm/ 4x8cm




Paper Tape

Width: 1.25cm/2.5cm/5cm/7.5cm/10cm



Self-adhesive bandage

Width: 2.5cm/5cm/7.5cm/10cm/15cm



Cotton ball

Gram weight: 0.3-4 g